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Instructions for use
MONDE | Orthodontics

Fundamentals

 Basic instructions for the use of MONDEAL Medical Systems GmbH (hereinafter referred to as MONDEAL) mini anchor screws and accessories in orthodontics. It is absolutely essential that all explanations and instructions in this manual are observed and fulfilled.

The implant may be inserted by only orthodontists, dentists, and oral surgeons as well as oral and maxillofacial surgeons.

If you are uncertain or have any questions, please contact us.

Information and symbols on labels

-  Item number
-  Quantity indicated on the packaging
-  LOT designation
-  Non-sterile
-  Sterilized using irradiation
-  Do not use if the packaging is damaged
-  Do not re-use
-  Manufacturer
-  Manufacturing date
-  Use before
-  Caution
-  Follow instructions for use
-  CE mark of conformity

Products

LOMAS – screws

Figure	Head shape	Diameter	Length	Colour
	Standard	1.5 mm	7–11 mm	Gold
		2.0 mm	7–11 mm	Grey
		2.3 mm	9 mm	Grey
	Hook	1.5 mm	7–11 mm	Gold
		2.0 mm	7–11 mm	Grey
		2.3 mm	9 mm	Grey
	Quad (horizontal, 0.018 × 0.025")	1.5 mm	7–11 mm	Gold
		2.0 mm	7–11 mm	Grey
		2.3 mm	9 mm	Grey
	Quad (horizontal, 0.022 × 0.028")	1.5 mm	7–11 mm	Green
		2.0 mm	7–11 mm	Blue
		2.3 mm	9 mm	Blue
	Quad (vertical, 0.018 × 0.025")	2 mm	7–11 mm	Grey
		2.3 mm	9 mm	Grey
	Quad (vertical, 0.022 × 0.028")	2 mm	7–11 mm	Blue
		2.3 mm	9 mm	Blue

LOMAS – VEGAS – Mini anchor screws

Figure	Diameter	Total length	Thread length	Colour
	2.0 mm	9 mm	7.0 mm	Grey
		11 mm	7.5 mm	
		13 mm	8.0 mm	
		15 mm	8.5 mm	
		17 mm	9.0 mm	
		20 mm	9.5 mm	
		25 mm	10.0 mm	
		30 mm	10.5 mm	
		35 mm	11.0 mm	
		40 mm	11.5 mm	

MONDEFIT screws

Figure	Diameter	Length	Colour
	1.5 mm	7–11 mm	Gold
	2.0 mm	7–11 mm	Grey
	2.3 mm	9 mm	Blue

Plates

Figure	Total length	Thickness	Colour
	29.5 mm	1.0 mm	Gold
	16.5 mm	1.0 mm	Grey

Plate screws

Figure	Diameter	Total length	Head shape	Colour
	2.0 mm	5 mm		Gold
		7 mm		Gold
	1.7 mm 2.0 mm	4 mm		Grey
		5 mm 7 mm		Grey

All MONDEAL implants are made of titanium. The screws are made of a titanium alloy Ti6Al4V ELI (ASTM F136, ISO 5832-3) and the plates of pure titanium (ASTM F67, ISO 5832-2). All titanium materials used are biocompatible, corrosion-resistant, and non-toxic in the biological environment. Titanium is completely non-magnetic.

Accessories

The LOMAS and MONDEFIT screws can be combined with commercially available orthodontic accessories such as rubbers/elastics, springs, and wires. With the MONDEFIT System, it is important to ensure that only specific MONDEAL MONDEFIT accessories – such as abutments, MONDEAL MONDEFIT plates, wires, fixing screws, hook locks and mobilisers – are used (see product catalogue).

The specific MONDEFIT accessories are made of special stainless implant steel (1.4441) in accordance with ISO 7153-1.

The impression cap used for taking impressions in the laboratory as well as the laboratory/manipulating implant are made of stainless steel 1.4305.

Incoming inspection

- Check the delivery for completeness and intactness immediately after receipt.
- Report any transport damage immediately.
- After delivery of the product, check the original packaging and the packaging seal for intactness. Packaging with a broken seal cannot be returned.
- Implant systems that are not packaged in a sterile condition must be removed from the packaging for reprocessing.
- Implants packed in sterile conditions can be applied directly from their sterile packaging without further pre-treatment.
- Sterile implants with (also inadvertently) opened or damaged sterile packaging as well as implants removed from packaging with an expired expiry date must be considered non-sterile and must be subjected to a validated reprocessing procedure before use (see chapter Notes on cleaning, disinfection and sterilisation).

Intended use

The implant is used for temporary anchorage in the bone to fix orthodontic appliances and aids for the treatment of malpositioned teeth.

Indications

LOMAS / MONDEFIT system:

- Distalisation and mesialisation of molars
- Molar anchorage, front anchorage
- Intrusion and extrusion
- Classification of retained teeth
- Alignment of molars
- Palatal expansion
- Temporary dentures for molars
- Change of the tooth position in the course of a pre-prosthetic treatment

Plates:

- Intrusion and extrusion of the teeth or groups of teeth
- Distalisation and mesialisation of molars
- Gap closing
- Three-dimensional possibility of tooth movements
- Movement of teeth or dental arch segments during skeletal anchorage
- Alignment of tilted molars

LOMAS Vegas screw:

- Change of the tooth position in the course of a pre-prosthetic treatment

- Molar anchorage, front anchorage
- Intrusion and extrusion
- Classification of retained teeth
- Alignment of molars

Clinical procedure

LOMAS / MONDEFIT screws

- Selection of the insertion site, screw length, and diameter using imaging techniques (measurement of bone and mucosa thickness).
- Local anaesthesia. Only surface anaesthesia if necessary.
- If necessary, pilot bore (12-12625 or 10-67517 depending on screw length). Mandibular and, as perforation of the compact, palatal. Manually with a suitable screwdriver or with motor system. 500-800 rpm under constant flushing with cooled, physiological saline solution.
- A separate drill should be used for each patient.
- Insertion with motor system: 10–20 N·cm recommended
- Insert the Mini Anchor Screw through the gingiva into the bone. Palatal: 90° to the occlusal plane possible, median and paramedian. Approx. 5–6 mm distance when using two mini anchor screws median. Paramedian: 3 mm from the midline on an imaginary transverse line between the palatal cusps of the first premolars. Orientation: third palatal fold.
- Immediate loading possible.
- The miniature anchor screw is removed in reverse order. No anaesthesia is necessary when unscrewing the screw.

MONDEFIT – plates

- The MONDEFIT system as a tandem screw (2 screws blocked with a MONDEPLATE) is used exclusively in the palate.
- The plates can be used both median and para-median.
- The MONDEFIT accessory is bent and attached after insertion of the MONDEFIT screw – either directly on the patient or in a laboratory using an impression cap (A33-VX1-0001) and a lab/manipulation implant (A33-VX1-0002) by means of an impression.
- During the bending process, make sure that the wire and the plate are held firmly and that the wire is not bent at the welding point. Otherwise, the welding point may break!

LOMAS Vegas – mini anchor bolts, short (9–13 mm) / medium (15–20 mm)

- Selection of the insertion site and screw length using imaging techniques and measurements
- Local anaesthesia.
- If necessary, pilot bore (10-67517 / A85-VA2-1500 - depending on AT) manually with a suitable screwdriver or with motor system. 500-800 rpm under constant flushing with cooled, physiological saline solution.
- Insertion with motor system: 10–20 N·cm recommended.
- Insert the Mini Anchor Screw through the gingiva into the bone with even rotational movements until the thread is screwed into the bone. Immediate loading possible.
- The mini anchor screw is removed by simply unscrewing it.

LOMAS Vegas – Mini anchor screws long (25–40 mm) – Insertion into the zygomatic

- Local anaesthesia.
- Vertical incision, upper jaw in the zygomatic of about 1 cm at about the height of the first molar.
- Subperiosteal dissection or preparation of the zygomatic bone.
- Drill with the 60 mm drill (A85-VA2-1500) as close as possible to the first molar buccal surface. The use of an angled hand piece facilitates drilling.
- Insert the screw with the corresponding blade manually (A33-18230) or with the angled hand piece (A33-18232 or A33-18237). A torque of 45 N·cm is recommended.
- Suture the incision – preferably with absorbable material. Immediately loadable.
- Removal is done by simply unscrewing the screw with a manual screwdriver or an angled hand piece.

Contraindications

- insufficient quantity or quality of bone mass/substance to prevent stable fixation with the device at or near the implant site
- lack of vascularisation
- poor oral hygiene
- other medical histories known to the user that have a negative influence on the healing process
- known allergies and/or hypersensitivity to implant materials
- general contraindications
- previous infections
- ulcers in the environment where the material is to be placed
- exposure to radioactive radiation or chemotherapy
- mental, physical, or neurological conditions that impair the patient's ability to heal postoperatively
- patients with insufficient ability and/or willingness to cooperate during the treatment phase

Possible complications

In most cases, possible complications are due to individual patient constitutions and general conditions rather than to the implant itself. These include:

- metal hypersensitivity or allergic reactions
- bone necrosis, osteoporosis
- insufficient revascularisation
- injuries to nerves and tendons
- irritation of the soft tissue
- superficial and deep infections
- implant loosening and implant loss
- pain and swelling

Taking into account the clinical condition and medical history of the patient, the treating physician must ensure that the use of implants is suitable for the individual case after weighing up a patient-specific risk/benefit assessment.

Mounting of sterile screws

- Remove the clip from the sterile packaging.
- Hold the clip on the grooved surface with two fingers, and pull the plastic cover to the side.
- Insert the insertion instrument into the clip and enclose the head of the screw. **Caution:** The head of the screw must sit securely in the receptacle of the insertion instrument.
- The screw can now be removed from the clip directly with the insertion instrument and inserted.

Warnings and precautions

- MONDEAL, as the manufacturer, recommends reading all available documents thoroughly before the first practical application and to contact users who already have practical experience with this type of treatment
- The products mentioned may only be used by qualified medical personnel with the appropriate training.
- The product must be handled and stored with care. Damage or scratches to the implant can significantly reduce the strength and fatigue resistance of the product.
- Products delivered non-sterile must be cleaned and sterilised by the user before use.
- All implant components are intended for single use and should never be reused.
- The patient must be advised not to manipulate the implant with the tongue in order to avoid loosening and loss of the implant. It is strongly recommended to insert the orthodontic appliance as soon as possible after the implants have been placed on the patient to avoid manipulation of the implant with the tongue and the resulting loosening of the implant.
- The patient must be instructed in postoperative hygiene.

MRI environment

The LOMAS, MONDEFIT and LOMAS VEGAS implants have **not** been tested for safety and compatibility in a computed tomography (CT) or magnetic resonance imaging (MRI) environment. For this reason, such use should be refrained from. The contraindicated use is the responsibility of the treating physician. The possible dangers include:

- Possible misinterpretation of the examination because of disturbing image artefacts.
- Heating or migration of the implant because of a reaction to the magnetic field.

Use of original products

Combination of different materials or systems is not permitted! Implant systems from MONDEAL must not be combined with systems from other manufacturers. If instruments are used that are not specifically designed for use with implants, the implant system may fail.

Special notes on disposable products

Depending on the product type, the first-time use of the product for its intended purpose may result in, among other things

- contamination of the product that can no longer be controlled safely
- material fatigue and material change
- non-visible damage, e.g. in the form of micro-cracks
- damage or scratches to the implant can significantly reduce the strength and fatigue resistance of the product.

Notes on suitable product selection

MONDEAL as a manufacturer does not recommend a certain surgical procedure for a specific patient. The operating surgeon is responsible for the choice of the appropriate implant in the case at hand.

The attending physician should be thoroughly familiar with the procedure through:

- careful review of the entire product documentation
- careful review of the current technical literature
- consultation with colleagues experienced in the use of this system
- practical exercises in the handling of the system, the surgical procedure, and the post-operative follow-up

In the case of complications, it may be necessary to remove the implant. Use the screwdriver provided for removal. Make sure that the screwdriver/screw head connection is exactly axial.

Reusability

The reuse of disposable products can compromise the structural integrity of the implants and increases the risk of functional failure. This can lead to injury to the patient. The reuse of disposable products also creates a risk of contamination (e.g. through germ transmission from patient to patient). This may result in injury and/or illness of the patient and/or user.

Instruments can be reused provided they are handled with care and are undamaged and unsoiled.

In case of misuse, MONDEAL excludes any liability.

Warranty

MONDEAL delivers only tested and faultless products to its customers. All of our products are designed and manufactured to meet the highest quality standards.

MONDEAL, as the manufacturer of the products, excludes any warranty claims and assumes no liability for direct or consequential damages resulting from:

- use for other purposes
- improper use, application, or handling
- lack of training
- improper preparation and sterilisation
- replacement of a MONDEAL product with a third-party product
- failure to observe the operating instructions

Notes on cleaning, disinfection, and sterilisation

Implants, instruments, containers, and accessories of the MONDEAL systems delivered in NON-STERILE condition must undergo suitable cleaning, disinfection, and sterilisation before use.

A description of the validated procedure for cleaning, disinfection, and sterilisation can be found in the **“Reprocessing instructions, cleaning, disinfection, and sterilisation”**.

Implants packed in **STERILE** condition can be applied directly from their sterile packaging without further pre-treatment.

Disposal

When disposing of packaging material and potentially infectious material (e.g. after metal removal), the national regulations and disposal guidelines must be observed.

Further information

Further information on the products (e.g. surgical techniques, care, cleaning, disinfection, and sterilisation) can be requested from your contact person. You can also find all information on the Internet at www.mondeal.de.

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