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Tel. +49 7463 99307 0
Fax: +49 7463 99307 33
E-Mail: mail@mondeal.de


Instruction for use MONDEAL Guiding and Comparative Instruments

Basic Information

The instruction for use for the Guiding and Comparative Instruments from MONDEAL Medical Systems GmbH (hereinafter referred to as MONDEAL).

It is essential that all the explanations and instructions in these instructions are observed and complied with.

Information and symbols on labels

	Item number
	Quantity
	Batch code
	Symbol for „UDI Code“
	Non-sterile
	Manufacturer
	Date of manufacture
	Caution
	Consult instructions for use
	CE- Conformity Declaration
	Symbol for „Prescription Only“ Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Notes Regarding the Delivered Goods

 The products may only be accepted when the manufacturer's labels and packaging are undamaged and unopened at the time of delivery.

The products are designed for repeated use. They are delivered **NON-STERILE** and must be appropriately cleaned, disinfected and sterilized before use. All packaging must be removed before cleaning

Product Material

The instruments are made of surgical steel according to ISO 7153-1 and ASTM F899 as well as Titanium according to ISO 5832-2 and ASTM F67. All materials are suitable for use as medical devices.

Variants

Guiding Instruments	
Item	Description
26-30201	Handle for drill guide only
26-30210	Drill guide for XS plates centric
26-30212	Drill guide for S-plates centric
26-30214	Drill guide for M / L plates centric
26-SG1-0001L	Boring jig for PRP left
26-SG1-0001R	Boring jig for PRP right
26-SG1-0002	Drilling guide for PRP plate
27-05142	Drill guide for Module 1 140 mm
27-05242	Drill guide for Module 2 140 mm
27-05342	Drill guide for Module 3 140 mm
27-05442	Drill guide for Module 4
34-10160	Boring bush without thread, with handle for BMR/MPS/Fracture-System

Comparative Instrumente	
Item	Description
10-60035	Depth comparison instrument Ø Screw 1.7 - 2.5 mm, Ø probe 1.2 mm Range 0 - 50 mm
26-30300	Depth comparison instrument Ø Screw 1.2 - 1.5 mm, Ø probe 0.8 mm, Range 0 - 30 mm

26-30311	Depth comparison instrument Ø Screw 1.7 - 2.5 mm, Ø probe 1.2 mm, Range 0 - 50 mm
27-06500	Depth comparison instrument for K-Wire 150 mm, M1 / M2 200 mm
27-06501	Depth comparison instrument for K-Wire 150 mm M3 200 mm
27-06502	Depth comparison instrument for K-Wire 150 200 mm M4 200 mm
27-06507	Screws and K-Wire comparison plate M1 / M2 / M3 150 mm
27-06508	Screws and K-Wire comparison plate M4 150 mm
27-08034	Depth comparison instrument Ø Screw 1.7 - 2.5 mm, Ø probe 1.2 mm, M1, Range 0 - 50 mm
27-08036	Depth comparison instrument Ø Screw 2.7 - 3.5 mm, Ø probe 1.7 mm, M2/M3, Range 0 - 80 mm
34-10155	Depth comparison instrument Ø Screw 1.7 - 2.5 mm, Ø probe 1.2 mm, Range 0 - 50 mm
27-06500	Depth comparison instrument for K-Wire 150 mm, M1 / M2 200 mm
27-06501	Depth comparison instrument for K-Wire 150 mm M3 200 mm
27-06502	Depth comparison instrument for K-Wire 200 mm M4 200 mm
84-SP0-0001	Screws comparison plate

Intended purpose

Guiding and comparative instruments are used to lead and compare shapes and geometries, vessels, tissues, organs, medical auxiliaries and implants. They are used in all surgical disciplines. They are reusable and can be used invasively in some cases

Indication

Guiding and comparative instruments are used for support in orthopedic surgery of the hand, foot, jaw, face and orthodontic surgery.

Contraindications

Guiding and comparative instruments are contraindicated in the following cases:

- Patients who do not meet the health requirements for surgical intervention.
- In patients in whom treatment with MONDEAL implants or prosthetic components is contraindicated.
- In patients who have an allergy or hypersensitivity to surgical steel or titanium. Information on contraindications specific to the implants can be found in the MONDEAL Instructions for Use for the respective systems

Possible Complications

In most cases, possible complications are caused due to patients' individual constitutions and general circumstances than due to the instruments themselves. These include, among others:

- Hypersensitivities to metals or other allergic reactions
- Injured nerves and tendons
- Irritation to soft tissue
- Superficial and deep infections
- Breakage of the instrument
- Pain and swelling

Taking into consideration the patient's clinical status and medical history, the attending physician must determine that the use of instruments is appropriate for the individual case in question, after weighing up the risk/benefit assessment for the specific patient.

Warnings and Precautions

- MONDEAL as a manufacturing company, recommends the user before the first practical application to thoroughly read all available documentation and contact users who have already had practical experience with this type of treatment procedure.
- The mentioned products must be used only by medical professionals with the appropriate training.
- The product must be handled and stored with care. Damage or scratches to the instrument may affect the safety of the product.
- Products damaged by transport, handling in the clinic or otherwise must not be used under any circumstances!
- Use the necessary care when using and storing the products.
- The products were developed and manufactured for a specific purpose. It must not be modified by the user other than as specified in these instructions for use.

Clinical procedure

Guiding Instruments

- Drilling aids and guides must be held firmly so that they do not slip when the drill is applied.

Comparative Instruments

- Before using the comparative instruments, check the zero position.
- The screws are held at the corresponding position/opening of the comparative plate using the screwdriver blade to determine the corresponding length.
- When checking the diameter, care must be taken to insert and remove the screws carefully to avoid jamming or stripping the screws.
- For depth comparators, determine the required screw length through the plate hole (i.e. with the plate in place).

Use of origin products

All system components have been designed and manufactured for a specific purpose and are therefore finely tuned to each other. No component may be modified by the user or replaced by any instrument or product from a third-party manufacturer, even if this resembles or even exactly matches the size or shape of the original product.

Storage

The instruments are to be stored dry and dust-protected at room temperature. There are no special requirements for the storage of the product.

Disposal

The products can be disposed of according to local requirements.

Return

Any return of products may only be sent back to us after disinfection/sterilization has been carried out and is clearly visible (appropriate packaging with sterile indicators, decontamination certificate, etc.).

The corresponding hygiene and business premises regulations must be observed.

Guarantee

MONDEAL delivers only tested and fault-free products to its customers. All our products are designed and manufactured to meet the highest quality standards. MONDEAL, as the manufacturer of the products, excludes all warranty claims, and assumes no liability for direct or consequential damages, resulting from:

- Unintended use.
- Improper use, application or handling
- Improper preparation and sterilization
- Not observing the instructions for use

Notes on cleaning, disinfection and sterilization

The products are delivered **NON-STERILE** and must be cleaned, disinfected and sterilized before use.

All packaging materials must be removed before reprocessing. Effective cleaning and disinfection are an indispensable prerequisite for effective sterilization.

As part of your responsibility for the sterility of the products, please ensure during application that, in principle, only sufficiently validated device-specific and product-specific methods for cleaning/disinfection and sterilization are used, that the equipment used (cleaning and disinfection unit, sterilizer) is regularly maintained and tested, and that the validated parameters are adhered to during every cycle.

Please also observe the applicable legal provisions in your country as well as the hygiene regulations of the doctor's practice or hospital. This applies in particular to the different specifications with regard to effective deactivation of prions.

Basic principles

If possible, a machine method should be used for cleaning and disinfection of the instruments (cleaning and disinfection unit/disinfector). Manual methods – even using an ultrasonic bath – should not be used due to the significantly lower effectiveness and possible damage. Pre-treatment must be carried out in both cases

Selection of cleaning agents, disinfectants and devices

When selecting the cleaning agents, disinfectants and devices used, it must be ensured at all stages that:

- These are suitable for the intended use (e.g., cleaning, disinfection, ultrasonic cleaning)
- The cleaning agents and disinfectants are non-protein-fixing (aldehyde-free)
- These are proven to be effective (e.g., VAH/DGHM or FDA approval or CE marking)
- The cleaning agents and disinfectants are suitable for the products and are compatible with the products

- The manufacturer's instructions, e.g., with regard to concentration, exposure time and temperature, are adhered to

When using **cleaning aids**, as well as during pre-cleaning, it must be ensured that:

- Only clean, lint-free cloths and/or soft brushes are used (never metal brushes or steel wool)
- If necessary, aids are used, such as cleaning pens, syringes, cannulas, bottle brushes for cannulated products or products with lumens.

For **drying** the products, MONDEAL recommends lint-free disposable cloths or medical compressed air.

With regard to water quality, MONDEAL recommends using demineralized and purified water (e.g. aqua purificata) for the cleaning, disinfection and rinsing steps.

Device for comparison depth:

During pre-cleaning, move the probe back and forth several times and disassemble if this is intended for cleaning.

Machine cleaning/disinfection

When selecting the cleaning agents and disinfectants, the chapter "Selecting cleaning agents, disinfectants and devices" must be observed. Likewise, the concentrations specified by the manufacturer of the cleaning agent and, if applicable, disinfectant must be adhered to, without fail.

When selecting the disinfectant, it must be ensured that the following phases are part of a cleaning process, in accordance with EN ISO 15883:

Phase	Temperature	Duration	Action
Rinsing	cold	1 min.	rinse with cold water
Cleaning	55°C	10 min.	add the cleaning agent
Neutralization	cold	2 min.	neutralize with cold water
Thermal disinfection	≥ 93°C	5 min.	with demineralized and purified water, do not add any additional cleaning agent
Drying	≥ 100°C	20 min.	drying process

When inserting the instruments into the cleaning and disinfection unit, make sure that the products do not touch each other and that they are aligned in such a way that no large liquid residues are able to remain on or in the product.

It is important to ensure that the products are well rinsed and that there is no foam residue.

Proof of the products' basic suitability for effective machine cleaning and disinfection was provided by an independent accredited testing laboratory using the "Miele PG 8535" cleaning and disinfection unit (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the cleaning agent "Neodisher mediclean forte" (Dr. Weigert GmbH & Co.

Inspection

Check all products for damage and contamination after cleaning and disinfecting and prior to sorting into containers and discard damaged and soiled products.

Packaging

Sort the cleaned and disinfected products into the sterilization containers and package them in disposable sterilization packs (single or double packaging) and/or sterilization containers that meet the following requirements:

- accord with DIN EN ISO/ANSI AAMI ISO 11607 and EN 868-2 to -10
- suitable for steam sterilization (temperature resistant up to min. 137 °C (279 ° F), sufficient vapor permeability)
- adequate protection of the products or sterilization packaging against mechanical damage
- regular maintenance according to the manufacturer's specifications (sterilization containers)

Sterilization

Use only the following listed sterilization methods for sterilization; other sterilization methods are not permitted.

All NON-STERILE products can be sterilized with steam in autoclaves (DIN EN 13060 or DIN EN 285).

Proof of the basic suitability of the instruments and implants for effective steam sterilization was provided by an independent accredited testing laboratory using the steam sterilizer "Lautenschläger ZentraCert" and the fractionated vacuum method. In the process, the method described above was observed:

Method	Fractionated or dynamic vacuum method	Flow, gravitation method
--------	---------------------------------------	--------------------------

IMPORTANT! READ CAREFULLY BEFORE USE; KEEP SAFE TO REFER TO AT A LATER POINT

Duration of exposure	≥ 4 min.	Not recommended
Temperature	132°C / 134°C	Not recommended
Drying time	> 20–30 Min.	Not recommended

MONDEAL recommends sterilization according to the fractionated or dynamic vacuum method. The use of the less effective gravitation method is only permissible if the fractionated vacuum process is unavailable; it also requires significantly longer sterilization times

If the user uses procedures other than those that are recommended, these must be validated by the user.

Ultimate responsibility for the validation of sterilization technology and equipment lies with the user.

Outside the US: The duration of sterilization can be extended to 18 minutes to comply with WHO recommendations and the Robert Koch Institute (RKI).

Do not use flash sterilization, hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization, or sterilization of thermolabile goods such as plasma or peroxide sterilization.

After sterilization, the products must be stored in the sterilization package dry and dust-free.

Note reporting obligations

All serious incidents that have occurred in connection with the product must be reported to the manufacturer, the distributor or the competent authority of the member state in which the user is established.

A report can be made at any time by sending an e-mail to qm@mondeal.de.

Additional information

You can request further information on the products (e.g., instructions for use, surgical techniques) from your contact person. In addition, you will find all the information on the Internet at www.mondeal.de.

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