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## Instruction for use for MONDEAL CMF-Accessories

### Basic information

The operating instructions refer to the CMF-Accessories from MONDEAL Medical Systems GmbH (hereinafter referred to as MONDEAL).

It is essential that all the explanations and instructions in these instructions are observed and complied with.

### Information and symbols on labels

	Item number
	Quantity
	Batch code
	Symbol for „UDI Code“
	Non-sterile
	Manufacturer
	Date of manufacture
	Caution
	Consult instructions for use
	CE- Conformity Declaration

### Notes Regarding the Delivered Goods

The products may only be accepted when the manufacturer's labels and packaging are undamaged and unopened at the time of delivery.

The products are intended for single use only and is not designed to be reused. It is delivered **NON-STERILE** and must be appropriately cleaned, disinfected and sterilized before use. All packaging must be removed before cleaning.

### Product Material

The products are made of pure titanium (ASTM F67, ISO 5832-2). This material is biocompatible, corrosion-resistant and non-toxic in a biological environment.

### Variants

#### Bending Template

68-00008	Template for BMR plate 34-40008 grey
68-00010	Template for BMR plate 34-40010 grey
68-00012	Template for BMR plate 34-40012 grey
68-00014	Template for BMR plate 34-40014 grey
68-00016	Template for BMR plate 34-40016 grey
68-00018	Template for BMR plate 34-40018 grey
68-00028	Template for BMR plate 34-40028 and 34-40048 grey
68-00032	Template for BMR plate 34-40032 and 34-40052 grey
68-00036	Template for BMR plate 34-40036 and 34-40056 grey
68-00040	Template for BMR plate 34-40040 and 34-40060 grey
68-00080	Template for BMR plate 34-40080 grey

#### Ringsplint

A30-07010	Maxillary Ringsplint 14-hole 0.6 mm grey
A30-07011	Maxillary Ringsplint 8-hole 0.6 mm grey
A30-07020	Mandibular Ringsplint 14-hole 0.6 mm grey
A30-07021	Mandibular Ringsplint 8-hole 0.6 mm grey

### Intended purpose

Bending Templates	Single-use template, for adapting all mandibular plates to the respective patient anatomy.
Ringsplint	Used for the treatment of traumatized or replanted teeth.

### Indication

Bending Templates	<ul style="list-style-type: none"> <li>Fractures,</li> <li>Osteotomies,</li> <li>Reconstructive procedures,</li> <li>Bone defects of the mandible.</li> </ul>
Ringsplint	The ring splints are used for: <ul style="list-style-type: none"> <li>stabilization after reduction or replantation (acute trauma cases)</li> <li>stabilization after intentional (planned) replantation</li> <li>stabilization in case of tooth luxation</li> </ul>

### Contraindications

Bending Templates	<ul style="list-style-type: none"> <li>General contraindications</li> <li>Known allergies and/or hypersensitivity.</li> </ul>
Ringsplint	<ul style="list-style-type: none"> <li>Replantation of deciduous teeth</li> <li>Known allergies and/or hypersensitivity</li> <li>Patients who are incapacitated and/or uncooperative during the treatment phase</li> </ul>

### Possible Complications

In most cases, possible complications are due to individual patient constitutions and general conditions rather than the implant/instruments themselves. Product related possible complications:

- Irritation of the soft tissue.
- Transient speech impairment.
- Use of excessive amounts of composite may cause hygiene problems.
- Inadequate composite fixation of the Ringsplint (e.g. due to insufficient drying/bonding) may result in loosening of the Ringsplint, which entails the risk of swallowing and aspiration of the device.
- Damage to the enamel during removal

Taking into account the patient's clinical condition and medical history, the treating physician must ensure that the use of implants is suitable for the individual case after weighing a patient-specific risk/benefit assessment.

### Warnings and Precautions

- MONDEAL as a manufacturing company, recommends the user before the first practical application to thoroughly read all available documentation and contact users who have already had practical experience with this type of treatment procedure.
- The mentioned products must be used only by medical professionals with the appropriate training.
- The product must be handled and stored with care. Damage or scratches to the instrument may affect the safety of the product.
- Products damaged by transport, handling in the clinic or otherwise must not be used under any circumstances!
- Use the necessary care when using and storing the products.
- The products were developed and manufactured for a specific purpose. It must not be modified by the user other than as specified in these instructions for use.
- Products are intended for single use on one patient. Products that have been used and removed once must be disposed of in accordance with local requirements. Reuse of products that have already been applied may compromise the structural integrity of the products and increases the risk of functional failure, which may result in patient injury. Reuse of devices that have already been used also creates a risk of contamination, e.g., through germ transmission from patient to patient. This may result in injury and/or illness to the patient and/or user. In case of disregard, MONDEAL excludes any liability.
- Products that have not come into direct contact with a patient may be reprocessed.
- Intraoral appliances must be secured against aspiration.

### MRT environment

The products have not been tested for safety and compatibility in a computer tomography (CT) or magnetic resonance imaging (MRI) environment. For this reason, such use is to be refrained from. The contraindicated use is the responsibility of the treating physician. The possible dangers include

- Possible misinterpretation of the examination due to disturbing image artefacts.
- Heating or migration of the products due to a reaction to the magnetic field.

## Instructions for Selecting the Appropriate Product

MONDEAL does not recommend specific procedures for specific patients or clinical conditions. The treating physician is solely responsible for choosing the appropriate treatment procedure and customizing an appropriate splint for the individual case. The follow-up treatment as well as the decision of when to remove the splint is the responsibility of the user. In general, the splints are designed for temporary fixation of traumatized/replanted teeth until the periodontium has sufficiently healed. Prolonged and rigid splinting may lead to adverse effects such as ankylosis or resorption. Splinting periods should be in accordance with clinical and radiological findings. The recommended splinting period is 2-4 weeks at most.

The treating dentist should become thoroughly familiarized with the procedure, for example by:

- Carefully studying all the product documentation
- Carefully reviewing the current scientific literature
- Consulting colleagues experienced in this field and with the use of this tooth splint
- Practicing handling of the products prior to its application

## Handling Instructions

### Bending Template

The bending templates are to be selected according to the size of the planned implant. The bending templates are held against the patient's freely prepared mandible and bent according to the patient's anatomical structure, using bending pliers and / or the fingers.

The finished bending templates serve as a template to bend the actual implant.

### Ringsplint

The adhesive surfaces of the Ringsplint are defined by the round openings of the splint. For sufficient stabilization of traumatized/replanted teeth, the Ringsplint should extend over at least two healthy teeth on both sides of the tooth/teeth to be stabilized. Be careful to position the Ringsplint with distance to the gingiva.

#### 1. Adapt the Ringsplint

Adapt the Ringsplint to match the patient's dentition and shape of the dental arch. Cut the splint to the desired length (e.g. by using a crown cutter or side cutter). Manually contour the splint to the curvature of the dental arch.

Caution: Do not repeatedly bend the splint in opposite directions since this may cause the splint to break.

#### 2. Clean the tooth surfaces

Thoroughly clean, rinse and air-dry the enamel surfaces.

#### 3. Etch the tooth surfaces

Etch the enamel surfaces thoroughly with ortho-phosphoric acid according to the instructions of the bonding agent manufacturer. Rinse with water and air-dry the tooth surfaces thoroughly after etching.

Caution: To allow appropriate bonded restoration, ensure to create an absolutely dry operating field.

#### 4. Bond the tooth surfaces

Apply the adhesive agent (bonding) according to the instructions of the bonding agent manufacturer. Light-cure the bonding agent as described by the manufacturer.

#### 5. Position the Ringsplint

Position the splint by hand or using tweezers. Ensure that the round openings of the Ringsplint match the etched/bonded sites. Take care to place the splint with distance to the gingiva.

#### 6. Fix the Ringsplint with flowable composite

Fix the Ringsplint to the teeth on both sides of the traumatized/replanted tooth/teeth by applying flowable composite to the rhomboid openings of the Ringsplint according to the instructions for use of the composite manufacturer. Light-cure the composite as described by the manufacturer. Reposition the traumatized/replanted tooth/teeth and bond it/them to the Ringsplint by applying flowable composite to the rhomboid openings of the Ringsplint according to the instructions of the composite manufacturer. Light-cure the composite as described by the manufacturer.

#### 7. Final trimming of the Ringsplint

Check the occlusion and check for sharp edges of the splint that could irritate the lips and soft tissues of the cheek. If necessary, modify the splint (recutting), especially at its ends.

Tip: To avoid/reduce possible soft tissue irritations, cover the ends of the splint with a small amount of composite.

Recommended splinting period: max. 2-4 weeks

## Removal

The decision of when to remove the splint is the responsibility of the treating physician. In general, the splint is designed for temporary fixation of traumatized/replanted teeth until sufficient periodontal healing has taken place. Prolonged and rigid splinting may lead to adverse effects such as ankylosis and resorption. Splinting periods should be in accordance with clinical and radiological findings of the individual case. The recommended splinting period is 2-4 weeks at most.

1. Remove the composite to the level of the splint using rotating instruments.
2. Hold the Ringsplint at one end and peel it off the tooth surfaces (manually or using tweezers).
3. Remove any residual composite (e.g. by using a dental curette) and refine the enamel surfaces with polishing disks. Be careful not to damage the enamel.
4. After removal of the splint, application of a fluoride-containing solution for remineralization of the etched enamel is highly recommended.

## Aftercare instruction

The treating physician has to instruct the patient on appropriate hygiene procedures according to standard dental practice. Depending on the clinical and radiological findings, the treating physician must advise the patient on possible restrictions regarding solid foods or sports following standard dental practice.

## Disposal

The products can be disposed of according to local requirements.

## Notes on cleaning, disinfection and sterilization

The products are delivered **NON-STERILE** and must be cleaned, disinfected and sterilized before use.

All packaging materials must be removed before reprocessing. Effective cleaning and disinfection are an indispensable prerequisite for effective sterilization.

As part of your responsibility for the sterility of the products, please ensure during application that, in principle, only sufficiently validated device-specific and product-specific methods for cleaning/disinfection and sterilization are used, that the equipment used (cleaning and disinfection unit, sterilizer) is regularly maintained and tested, and that the validated parameters are adhered to during every cycle.

Please also observe the applicable legal provisions in your country as well as the hygiene regulations of the doctor's practice or hospital. This applies in particular to the different specifications with regard to effective deactivation of prions.

## Basic principles

If possible, a machine method should be used for cleaning and disinfection of the instruments (cleaning and disinfection unit/disinfector). Manual methods – even using an ultrasonic bath – should not be used due to the significantly lower effectiveness and possible damage. Pre-treatment must be carried out in both cases

## Selection of cleaning agents, disinfectants and devices

When selecting the cleaning agents, disinfectants and devices used, it must be ensured at all stages that:

- These are suitable for the intended use (e.g., cleaning, disinfection, ultrasonic cleaning)
- The cleaning agents and disinfectants are non-protein-fixing (aldehyde-free)
- These are proven to be effective (e.g., VAH/DGHM or FDA approval or CE marking)
- The cleaning agents and disinfectants are suitable for the products and are compatible with the products
- The manufacturer's instructions, e.g., with regard to concentration, exposure time and temperature, are adhered to

When using **cleaning aids**, as well as during pre-cleaning, it must be ensured that:

- Only clean, lint-free cloths and/or soft brushes are used (never metal brushes or steel wool)
- If necessary, aids are used, such as cleaning pens, syringes, cannulas, bottle brushes for cannulated products or products with lumens.

For **drying** the products, MONDEAL recommends lint-free disposable cloths or medical compressed air.

With regard to water quality, MONDEAL recommends using demineralized and purified water (e.g. aqua purificata) for the cleaning, disinfection and rinsing steps.

## Machine cleaning/disinfection

When selecting the cleaning agents and disinfectants, the chapter "Selecting cleaning agents, disinfectants and devices" must be observed. Likewise, the concentrations specified by the manufacturer of the cleaning agent and, if applicable, disinfectant must be adhered to, without fail.

When selecting the disinfector, it must be ensured that the following phases are part of a cleaning process, in accordance with EN ISO 15883:

## IMPORTANT! READ CAREFULLY BEFORE USE; KEEP SAFE TO REFER TO AT A LATER POINT

Phase	Temperature	Duration	Action
Rinsing	cold	1 min.	rinse with cold water
Cleaning	55°C	10 min.	add the cleaning agent
Neutralization	cold	2 min.	neutralize with cold water
Thermal disinfection	≥ 93°C	5 min.	with demineralized and purified water, do not add any additional cleaning agent
Drying	≥ 100°C	20 min.	drying process

When inserting the implants and instruments into the cleaning and disinfection unit, make sure that the products do not touch each other and that they are aligned in such a way that no large liquid residues are able to remain on or in the product.

It is important to ensure that the products are well rinsed and that there is no foam residue.

Proof of the products' basic suitability for effective machine cleaning and disinfection was provided by an independent accredited testing laboratory using the "Miele PG 8535" cleaning and disinfection unit (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the cleaning agent "Neodisher mediclean forte" (Dr. Weigert GmbH & Co.

### Inspection

Check all products for damage and contamination after cleaning and disinfecting and prior to sorting into containers and discard damaged and soiled products.

### Packaging

Sort the cleaned and disinfected products into the sterilization containers and package them in disposable sterilization packs (single or double packaging) and/or sterilization containers that meet the following requirements:

- accord with DIN EN ISO/ANSI AAMI ISO 11607 and EN 868-2 to -10
- suitable for steam sterilization (temperature resistant up to min. 137 °C (279 ° F), sufficient vapor permeability)
- adequate protection of the products or sterilization packaging against mechanical damage
- regular maintenance according to the manufacturer's specifications (sterilization containers)

### Sterilization

Use only the following listed sterilization methods for sterilization; other sterilization methods are not permitted.

All NON-STERILE products can be sterilized with steam in autoclaves (DIN EN 13060 or DIN EN 285).

Proof of the basic suitability of the instruments and implants for effective steam sterilization was provided by an independent accredited testing laboratory using the steam sterilizer "Lautenschläger ZentraCert" and the fractionated vacuum method. In the process, the method described above was observed:

Method	Fractionated or dynamic vacuum method	Flow, gravitation method
Duration of exposure	≥ 4 min.	Not recommended
Temperature	132°C / 134°C	Not recommended
Drying time	> 20–30 Min.	Not recommended

MONDEAL recommends sterilization according to the fractionated or dynamic vacuum method. The use of the less effective gravitation method is only permissible if the fractionated vacuum process is unavailable; it also requires significantly longer sterilization times

If the user uses procedures other than those that are recommended, these must be validated by the user.

Ultimate responsibility for the validation of sterilization technology and equipment lies with the user.

Outside the US: The duration of sterilization can be extended to 18 minutes to comply with WHO recommendations and the Robert Koch Institute (RKI).

Do not use flash sterilization, hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization, or sterilization of thermolabile goods such as plasma or peroxide sterilization.

After sterilization, the products must be stored in the sterilization package dry and dust-free.

### Note reporting obligations

All serious incidents that have occurred in connection with the product must be reported to the manufacturer, the distributor or the competent authority of the member state in which the user is established.

A report can be made at any time by sending an e-mail to [qm@mondeal.de](mailto:qm@mondeal.de).

### Additional information

You can request further information on the products (e.g., instructions for use, surgical techniques) from your contact person. In addition, you will find all the information on the Internet at [www.mondeal.de](http://www.mondeal.de).

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